

Hamilton Dilutors/Dipsenser, syringe- and valve-drives have a serial interface to control the devices. The Software ML-Script is made to simplify controlling the instruments.

You can visually create your own methods with a windows software. The program generates the requested commands, send it to the devices and control their response and error status. If the device should run wo/PC you can download the method into the device's memory.

Supported instruments

- ML500-Series (1- or 2-syringe dispenser or dilutors)
- PSD/2-syringe drives
- PSD/3-syringe drives

{xtypo_download}Datasheet ([View](#) / [Download](#)){/xtypo_download}

{xtypo_download}Manual ([View](#) / [Download](#)) {/xtypo_download}